Cat's Pur



Uses

- Used in DCR/DCT to retract skin, muscle
- Dacryocystectomy—DCT—entire lacrimal sac is removed.

Contraindications for DCR

- 1. In elderly due to atrophied sac
- 2. In children less than 4 years as the bone growth is less
- 3. Infections of the sac like TB, syphilis
- 4. Malignancy of the sac
- 5. Uncorrected nasal pathology like DNS, nasal polyp
- 6. Atrophic rhinitis.

Lacrimal Syringing Cannula



Uses

Used to check patency of lacrimal drainage system.

Interpretation of syringing

- 1. While syringing if the fluid comes through the same punctum, it suggests canalicular block.
- 2. While syringing if the fluid comes from the opposite punctum, it suggests nasolacrimal duct or common canalicular block.

Citelli's Bone Punch



Uses

- To remove the bone in external DCR

 Bones removed are frontal process of the maxilla, lacrimal bone
- Size of the opening should be at least 10 mm.

Cryo Pencil



Uses

- For ICCE
- Filled with freon.

Indications for ICCE

Grossly subluxated/dislocated lens

Other uses of cryo in ophthalmology

Distichiasis, cyclocryotherapy for absolute glaucoma, to seal the retinal holes in retinal detachement surgery, treatment of retinoblastoma.

Joule Thomson effect

When a non-ideal gas suddenly expands from a high pressure to low pressure, there is a drop in the temperature. This principle is used in cryoapplication.

Manipal Manual of Instruments		

Manipal Manual of Instruments

Contents

Foreword	vii
Preface	ix
Contributors	xi
Section 1	
Common Instruments	1
Section 2	
ENT Instruments	13
Section 3	
Ophthalmic Instruments	67
Section 4	
General Surgical Instruments	95
Section 5	
Obstetrical and	
Gynaecological Instruments	141
Section 6	
Orthopaedic Instruments,	
Implants and Prostheses	187
Section 7	
Anaesthesia Instruments	221
Index	277